



## **API Guide**

Contents

*REST Client Authentication Mechanism ..... 1*

*Ordr SCE REST APIs ..... 1*

*REST APIs - Sample Queries and Output ..... 4*

## REST Client Authentication Mechanism

REST server supports three modes of authentication, as follows:

- Authentication through an AD server
- Authentication through LDAP server
- Authentication using credentials configured on the Ordr SCE

The AD/LDAP authentication servers' configuration page is located at **System > Service Integration > Internal Services > External Authentication** in Ordr SCE user interface.

The REST server authentication mode configuration page is located at **System > Service Integration > Internal Services > SCE API** in Ordr SCE user interface.

### Note:

The response to a request may not contain all the data, because the REST server supports paginated response. In paginated response, the response includes metadata that contains the next-link the REST client must submit in order to fetch successive data for the request it had originally initiated.

## Ordr SCE REST APIs

HTTP Method	Uniform Resource Identifier (URI)	API Description	Output Format
GET	/Rest/Devices	<a href="#">Fetches information of all devices in the system.</a>	JSON
GET	/Rest/Devices/version	<a href="#">Retrieves the current version of the Devices API supported by Ordr SCE.</a>	Plain text
	/Rest/Devices?limit=<value>	Controls the count of devices in response. For example: Devices?limit=10	
GET	/Rest/Devices?vulnIds	<a href="#">Fetches the list of devices having specific vulnerabilities.</a> <b>Note:</b> Values can be comma separated list. For example: /Devices?vulnIds==FDA-173239,CVE-2020-0601	JSON
	/Rest/Devices?iot=true	Fetches all IoT devices.	
	/Rest/Devices?non_iot=true	Fetches all Non-IoT devices.	
	/Rest/Devices?openPorts=true	Includes open ports in devices information response.	

	/Rest/Devices?weakPassword=true	Includes open ports with weak password in devices information response.	
	/Rest/Devices?tenantGuid=<tenantGuid>	Uses APIs for a particular tenant.	
	/Rest//Devices?include=subcategory	Includes optional information in response. <b>Note:</b> Values can be comma separated list.	
GET	/Rest/Devices?mac=<mac-address>	<a href="#">Fetches the device information for the given MAC address.</a>  The software's information (third-party apps, OS patches, and anti-virus products) is displayed in response if the device is Windows and WinRM enabled.	JSON
GET	/Rest/Devices?ip=<ip-address>	<a href="#">Fetches the device information for the given IP address.</a>  The software's information (third-party apps, OS patches, and anti-virus products) is displayed in response if the device is Windows and WinRM enabled.	JSON
GET	/Rest/Devices?group=<group-name>	<a href="#">Fetches all the devices belonging to the given group. For example, medical devices group.</a>	JSON
GET	/Rest/Devices?appName=<application-name>	<a href="#">Fetches all devices talking to an application.</a>	JSON
GET	/Rest/Devices?riskState=<state>	<a href="#">Fetches all the devices whose risk state matches the specified state.</a>  Possible values for risk state are Critical, High, Medium, Low, and Normal.	JSON
GET	/Rest/Devices?connStatus=ONLINE_IN_LAST_24_HRS	<a href="#">Fetches new devices that showed up in the last 24 hours.</a>	JSON
GET	/Rest/Devices?connStatus=ONLINE_IN_LAST_WEEK	<a href="#">Fetches new devices that showed up within last week.</a>	JSON
GET	/Rest/Devices?connStatus=ONLINE&startTime=1569737345520&endTime=1569823745520	<a href="#">Fetches new devices that showed up between start time and end time.</a>  <b>Note:</b> You should specify start time and end time in milliseconds (since Jan 1, 1970, UTC).	JSON
GET	/Rest/Devices?connStatus=ONLINE	<a href="#">Fetches devices that are currently online.</a>	JSON
GET	/Rest/Devices?connStatus=OFFLINE_IN_LAST_24_HRS	<a href="#">Fetches devices that went offline in the last 24 hours.</a>	JSON
GET	/Rest/Devices?connStatus=OFFLINE_IN_LAST_WEEK	<a href="#">Fetches devices that went offline within last week.</a>	JSON
GET	/Rest/Devices?connStatus=OFFLINE&startTime=	<a href="#">Fetches devices that went offline between start time and end time.</a>	JSON

	1569737345520&endTime=1569823745520	<b>Note:</b> You should specify start time and end time in milliseconds (since Jan 1, 1970, UTC).	
	/Rest/DeviceUtil?mac=<mac address>&startTime=<unix epoch timestamp>&endTime=<unix epoch timestamp>	<a href="#">Fetches utilization information for a given device between start time and end time.</a>	
	/Rest/DeviceUtil?deviceType=<type>&startTime=<unix epoch timestamp>&endTime=<unix epoch timestamp>	<a href="#">Fetches utilization information for all devices of a given type (for example: patient monitor, infusion pump, or ultrasound/mri) between start time and end time.</a>	
	/Rest/DeviceUtil?profileName=<order profile name>&startTime=<unix epoch timestamp>&endTime=<unix epoch timestamp>	<a href="#">Fetches utilization information for all devices of a given profile between start time and end time.</a>	
GET	/Rest/Devices?connStatus=OFFLINE	<a href="#">Fetches devices that are currently offline.</a>	JSON
GET	/Rest/Devices?sensorName=test	<a href="#">Fetches devices that are currently under the purview of the test sensor.</a>	JSON
GET	/Rest/Devices?sensorIp=192.168.101.1	<a href="#">Fetches devices that are currently under the purview of the sensor with 192.168.101.1 IP.</a>	JSON
GET	/Rest/Devices?diskEncrypted=<status>	<a href="#">Fetches all devices which has its disk encryption status matches given status.</a> Status can be 'true' or 'false' (currently Windows-only devices using WinRM gathered data. Also return devices only if disk encryption status present in database).	JSON
GET	/Rest/Devices?biosPassword=true	<a href="#">Fetches all devices which has its BIOS password status matches given status.</a> Status can be 'true' or 'false' (currently Windows-only devices using WinRM gathered data. Also return devices only if BIOS password status present in database).	JSON
GET	/Rest/Devices?softwareInstalled=<softwareName>	<a href="#">Fetches all devices which has the specified software installed.</a> This will do a substring comparison of the passed argument. Example software name is 'Mozilla' (currently Windows-only devices using WinRM gathered data. Returns devices only if its software information present in database).	JSON
GET	/Rest/Devices?patchInstalled=<HotfixId>	<a href="#">Fetches all devices which has the specified hotfix installed.</a>	JSON

		This will do full string comparison of the passed argument. Example hotfix ID is KB4534132 (currently Windows-only devices using WinRM gathered data. Returns devices only if patch info present in database).	
GET	/Rest/Flows/version	<a href="#">Retrieves the current version of the Flows API supported by Ordr SCE.</a>	Plain text
GET	/Rest/Flows?srcIp=<ip-address>	<a href="#">Fetches all the flows whose source is the specified IP address.</a>	JSON
GET	/Rest/Flows?dstIp=<ip-address>	<a href="#">Fetches all the flows whose destination is the specified IP address.</a>	JSON
GET	/Rest/Flows?srcMac=<mac-address>	<a href="#">Fetches all the flows whose source is the specified device.</a>	JSON
GET	/Rest/Flows?dstMac=<mac-address>	<a href="#">Fetches all the flows whose destination is the specified device.</a>	JSON
	/Rest/Flows?limit=<value>	Controls the count of flows in response. For example: /Flows?limit=10	
GET	/Rest/Applications/version	<a href="#">Retrieves the current version of the Applications API supported by Ordr SCE.</a>	Plain text
GET	/Rest/Applications?mac=<mac-address>	<a href="#">Fetches all the applications used by the device for the given device MAC.</a>	JSON
GET	/Rest/Applications?ip=<ip-address>	<a href="#">Fetches all the applications used by the device for the given device IP.</a>	JSON
GET	/Rest/SecurityAlarms	<a href="#">Fetches all the security incidents detected by the system.</a>	JSON
GET	/Rest/SecurityAlarms?mac=<mac-address>	<a href="#">Fetches all the security incidents for the device for the given device MAC.</a>	JSON
GET	/Rest/SecurityAlarms?ip=<ip-address>	<a href="#">Fetches all the security incidents for the device for the given device IP.</a>	JSON
GET	/Rest/SecurityAlarms?category=<category-type>	<a href="#">Fetches all the security incidents for the given alarm category.</a>	JSON
GET	/Rest/SecurityAlarms/Summary	<a href="#">Provides a summary of the security alarms by category.</a>	JSON
	/Rest/SecurityAlarms?limit=<value>	Controls the count of alarms in response. For example: /SecurityAlarms?limit=10	
GET	/Rest/NetworkDevices	<a href="#">Fetches all network equipments.</a>	JSON

## REST APIs - Sample Queries and Output

**Description** : Retrieve the current version of the Devices API supported by Ordr SCE.

**Request URI** : https://<Ordr\_SCE>/Rest/Devices/version

**Output** : Current supported version of /Rest/Devices API is 1.0.

**Description** : Fetch device information for the given MAC address.

**Request URI** : <https://192.168.104.182/Rest/Devices?mac=00:50:56:6A:42:46>

**Output** :

```
{
  "MetaData":{
    "Count":1
  },
  "Devices":[
    {
      "MacAddress":"<mac-address>",
      "IpAddress":"<IP-address>",
      "Group":"Medical Devices",
      "Profile":"GE-LOGIQ700-Ultrasound",
      "MfgName":"GEMedica",
      "LongMfgName":"G.E. Medical Systems",
      "Vlan":204,
      "ModelNameNo":"LOGIQ 700",
      "RiskState":"NORMAL",
      "DeviceType":"Ultrasound",
      "SerialNo":"<serial-number>",
      "DeviceDescr":"Ultrasound",
      "Subnet":"10.200.204.0/24",
      "SwVersion":"R6.1"
      "softwareInfo": {
        "ThirdPartyApps": [{
          "Name": "Mozilla Firefox 65.0.1 (x64 en-US)",
          "Version": "65.0.1",
          "Vendor": "Mozilla",
          "InstallDate": null
        }],
        "OsPatches": [{
          "HotfixId": "KB4534132",
          "InstalledOn": "1581580800000",
          "Description": "Update"
        }],
        "AvProducts": [{
          "displayName": "Windows Defender",
          "ProtectionState": "ACTIVE",
          "IsUpToDate": "true",
```

```

        "UpdateTime": "true",
        "pathToSignedProductExe": "windowsdefender://"
    }
}
}

```

**Description** : Fetch the device information for the given IP address.

**Request URI** : [https://<Ordr\\_SCE>/Rest/Devices?ip=192.168.53.4](https://<Ordr_SCE>/Rest/Devices?ip=192.168.53.4)

**Output** : [Sample output](#)

**Description** : Fetch information of all devices in the system.

**Request URI** : [https://<Ordr\\_SCE>/Rest/Devices](https://<Ordr_SCE>/Rest/Devices)

**Output** :

```

{
  "MetaData": {
    "Count": 100,
    "next": "/Rest/Devices?clientMacToken=-7079632916954617042"
  },
  "Devices": [
    {
      "MacAddress": "<mac>",
      "IpAddress": "<ip>",
      "Group": "Network Devices",
      "Profile": "Cisco-WS-C3560X-24T-Catalyst Switch",
      "MfgName": "Cisco",
      "LongMfgName": "Cisco Systems, Inc",
      "Vlan": 2,
      "ModelNameNo": "WS-C3560X-24T",
      "RiskState": "NORMAL",
      "DeviceType": "Catalyst Switch",
      "SerialNo": "<serial>",
      "DeviceDescr": "Catalyst Switch",
      "SwVersion": "12.2(55)SE10",
      "OsVersion": "C3560E Software",
      "OsType": "Cisco IOS",
      "endpointType": "NONIOT_ENDPOINT",
    }
  ]
}

```



```
"knownVulnRiskState": "NORMAL",
"noOfPorts": 43,
"ports": [
  {
    "name": "GigabitEthernet0/12",
    "hardware": "Gigabit Ethernet",
    "type": "TRUNK",
    "vlan": "1",
    "remoteNwEquipIp": "10.200.201.8",
    "remoteNwEquipMac": "<mac>",
    "remoteNwEquipPort": "FastEthernet1/0/24",
    "remoteNwEquipName": "cisco_mgmt.not",
    "remoteNwEquipManufacturer": "Cisco Systems, Inc",
    "remoteNwEquipModelNo": "WS-C3750-24PS-S",
    "remoteNwEquipSwVersion": "12.2(55)SE7"
  }
],
"alarmCount": 0,
"riskScore": 0,
"firstSeen": "2020-02-19 05:54:40 GMT",
"lastSeen": "2020-03-15 06:03:10 GMT",
"classificationState": "Classified",
"sensorName": "reports-dpvm-ss48",
"sensorIp": "172.18.10.16",
"connStatus": "ONLINE"
},
{
  "MacAddress": "<mac>",
  "Group": "Medical Devices",
  "Profile": "Philips-Patient Monitoring",
  "MfgName": "PhilipsP",
  "LongMfgName": "Philips Patient Monitoring",
  "Vlan": 777,
  "ModelNameNo": "",
  "RiskState": "NORMAL",
  "DeviceType": "Patient Monitoring",
  "DeviceDescr": "Patient Monitoring",
  "OsType": "Linux Embedded RTOS",
  "endpointType": "IOT_ENDPOINT",
```

```

    "knownVulnRiskState": "NORMAL",
    "accessType": "WIRED",
    "nwEquipInterface": "52",
    "nwEquipHostname": "Aruba-2930F-48G-4SFPP",
    "nwEquipScrapeIp": "10.200.201.39",
    "alarmCount": 0,
    "riskScore": 0,
    "firstSeen": "2020-02-28 23:02:55 GMT",
    "lastSeen": "2020-02-28 23:02:55 GMT",
    "classificationState": "Classified",
    "sensorName": "reports-dpvm-ss48",
    "sensorIp": "172.18.10.16",
    "connStatus": "OFFLINE"
  },
  <<<<<SNIP>>>>>,
  {
    "MacAddress": "<mac>",
    "Group": "Mobile Phones and Tablets",
    "Profile": "Samsung-Galaxy Note9-Phone",
    "MfgName": "Samsung",
    "LongMfgName": "Samsung",
    "Vlan": 2,
    "ModelNameNo": "Galaxy Note9",
    "RiskState": "NORMAL",
    "DeviceType": "Phone",
    "DeviceDescr": "Phone",
    "OsType": "Android",
    "fqdn": "Galaxy-Note9.hq.ordr.net",
    "dhcpHostname": "Galaxy-Note9",
    "endpointType": "NONIOT_ENDPOINT",
    "knownVulnRiskState": "NORMAL",
    "alarmCount": 0,
    "riskScore": 0,
    "firstSeen": "2020-02-20 17:38:16 GMT",
    "lastSeen": "2020-03-06 01:59:33 GMT",
    "classificationState": "Classified",
    "sensorName": "dc13-dpvm-ss72",
    "sensorIp": "192.168.104.86",
    "connStatus": "OFFLINE"
  }

```

```

    }
  ]
}

```

**Description** : Fetch the list of devices having specific vulnerabilities.

**Request URI** : `https://<SCE-IP>/Rest/Devices?vulnIds=FDA-173239,CVE-2020-0601`

**Output** : [Sample output](#)

**Description** : Fetch information of all devices talking to an application.

**Request URI** : `https://<Ordr_SCE>/Rest/Devices?appName=udp.bacnet`

**Output** : [Sample output](#)

**Description** : Fetch devices that are currently online.

**Request URI** : `https://<sce ip>/Rest/Devices?connStatus=ONLINE`

**Output** : [Sample output](#)

**Description** : Fetch utilization information for a given device between start time and end time.

**Request URI** : `https://<sce ip>/Rest/DeviceUtil?mac=<mac address>&startTime=<unix epoch timestamp>&endTime=<unix epoch timestamp>`

**Output** :

```

{
  "MetaData":{"Count":1},
  "DeviceUtilRecords":
  [
    {
      "MacAddress":"48:0F:CF:48:40:BD",
      "MfgName":"GE MEDICAL SYSTEMS",
      "DeviceType":"MRI",
      "ModelNameNo":"Signa Pioneer",
      "UtilPercent":16
    }
  ]
}

```

**Description** : Fetch utilization information for all devices of a given type (for example: patient monitor, infusion pump, or ultrasound/mri) between start time and end time.

**Request URI** : `https://<sce ip>/Rest/DeviceUtil?deviceType=<type>&startTime=<unix epoch timestamp>&endTime=<unix epoch timestamp>`

**Output** :

```
{
  "MetaData":{"Count":4},
  "DeviceUtilRecords":
  [
    {
      "MacAddress":"00:80:17:3D:D6:40",
      "MfgName":"Hitachi Medical Corporation",
      "DeviceType":"MRI",
      "ModelNameNo":"Oasis",
      "UtilPercent":40
    },
    {
      "MacAddress":"00:1B:21:02:1C:46",
      "MfgName":"GE MEDICAL SYSTEMS",
      "DeviceType":"MRI",
      "ModelNameNo":"Signa HDxt",
      "UtilPercent":44
    },
    {
      "MacAddress":"C8:D3:FF:BA:7D:70",
      "MfgName":"Philips Medical Systems",
      "DeviceType":"MRI",
      "ModelNameNo":"Achieva",
      "UtilPercent":36
    },
    {
      "MacAddress":"00:0E:0C:F5:FB:1C",
      "MfgName":"GE MEDICAL SYSTEMS",
      "DeviceType":"MRI",
      "ModelNameNo":"Signa HDxt",
      "UtilPercent":18
    }
  ]
}
```

```
}
```

**Description** : Fetch utilization information for all devices of a given profile between start time and end time.

**Request URI** : `https://<sce ip>/Rest/DeviceUtil?profileName=<ordr profile name>&startTime=<unix epoch timestamp>&endTime=<unix epoch timestamp>`

**Output** :

```
{
  "MetaData":{"Count":2},
  "DeviceUtilRecords":
  [
    {
      "MacAddress":"48:0F:CF:48:40:BD",
      "MfgName":"GE MEDICAL SYSTEMS",
      "DeviceType":"MRI",
      "ModelNameNo":"Signa Pioneer",
      "UtilPercent":16
    },
    {
      "MacAddress":"30:9C:23:41:9E:FE",
      "MfgName":"GE MEDICAL SYSTEMS",
      "DeviceType":"MRI",
      "ModelNameNo":"Signa Pioneer",
      "UtilPercent":0
    }
  ]
}
```

**Description** : Fetch devices that are currently offline.

**Request URI** : `https://<sce ip>/Rest/Devices?connStatus=OFFLINE`

**Output** : [Sample output](#)

**Description** : Fetch new devices that showed up between start time and end time.

**Request URI** : `https://<sce ip>/Rest/Devices?connStatus=ONLINE&startTime=1569737345520&endTime=1569823745520`

**Output** : [Sample output](#)

**Description** : Fetch new devices that showed up in the last 24 hours.

**Request URI** : `https://<sce ip>/Rest/Devices?connStatus=ONLINE_IN_LAST_24_HRS`

**Output** : [Sample output](#)

**Description** : Fetch new devices that showed up within last week.

**Request URI** : `https://<sce ip>/Rest/Devices?connStatus=ONLINE_IN_LAST_WEEK`

**Output** : [Sample output](#)

**Description** : Fetch devices that went offline in the last 24 hours.

**Request URI** : `https://<sce ip>/Rest/Devices?connStatus=OFFLINE_IN_LAST_24_HRS`

**Output** : [Sample output](#)

**Description** : Fetch devices that went offline within last week.

**Request URI** : `https://<sce ip>/Rest/Devices?connStatus=OFFLINE_IN_LAST_WEEK`

**Output** : [Sample output](#)

**Description** : Fetch all the devices whose risk state matches the specified state. Possible values for risk state are Critical, High, Medium, Low, and Normal.

**Request URI** : `https://<Ordr_SCE>/Rest/Devices?riskState=MEDIUM`

**Output** : [Sample output](#)

**Description** : Fetch devices that are currently under the purview of the test sensor.

**Request URI** : `https://<sce ip>/Rest/Devices?sensorName=test`

**Output** : [Sample output](#)

**Description** : Fetch devices that are currently under the purview of the sensor with 192.168.101.1 IP.

**Request URI** : `https://<sce ip>/Rest/Devices?sensorIp=192.168.101.1`

**Output** : [Sample output](#)

**Description** : Fetch all devices which has its disk encryption status matches given status (true or false).

**Request URI** : `https://<Ordr_SCE>/Rest/Devices?diskEncrypted=<status>`

**Output** : [Sample output](#)

**Description** : Fetch all devices which has its BIOS password status matches given status (true or false).

**Request URI** : `https://<Ordr_SCE>/Rest/Devices?biosPassword=true`

**Output** : [Sample output](#)

**Description** : Fetch all devices which has the specified software installed.

**Request URI** : `https://<Ordr_SCE>/Rest/Devices?softwareInstalled=<softwareName>`

**Output** : [Sample output](#)

**Description** : Fetch all devices which has the specified hotfix installed.

**Request URI** : `https://<Ordr_SCE>/Rest/Devices?patchInstalled=<HotfixId>`

**Output** : [Sample output](#)

**Description** : Fetch all the devices belonging to the given group. For example, medical devices group.

**Request URI** : `https://<Ordr_SCE>/Rest/Devices?group=Industrial Devices`

**Output** : [Sample output](#)

**Description** : Retrieve the current version of the Flows API supported by Ordr SCE.

**Request URI** : `https://<Ordr_SCE>/Rest/Flows/version?srcIP`

**Output** : Current supported version of /Rest/Flows API is 1.0

**Description** : Fetch all the flows whose source is the specified IP Address.

**Request URI** : `https://<Ordr_SCE>/Rest/Flows?srcIp=10.22.22.176`

**Output** :

```
{
  "MetaData": {
    "Count": 2
  },
  "Flows": [
```

```
{
  "vectorGuid": "a-13913-1700256025",
  "behaviorState": "NORMAL",
  "srcIp": "10.200.204.9",
  "dstIp": "192.168.101.145",
  "srcPort": 49866,
  "dstPort": 104,
  "ipProto": 6,
  "rxBytes": 0,
  "txBytes": 60,
  "rxPkts": 0,
  "txPkts": 1,
  "lastSeenTimestamp": 1537417950408,
  "external": false,
  "appName": "https",
  "remoteProfile": "Local-IP-Profile",
  "alarms": [{
    "alarmHash": "<alarmHash>",
    "category": "<category>",
    "categoryType": "<categoryType>"
  }]
},
{
  "vectorGuid": "a-13917-1700256025",
  "behaviorState": "NORMAL",
  "srcIp": "10.200.204.9",
  "dstIp": "192.168.101.241",
  "srcPort": 60496,
  "dstPort": 53,
  "ipProto": 17,
  "rxBytes": 73,
  "txBytes": 57,
  "rxPkts": 1,
  "txPkts": 1,
  "lastSeenTimestamp": 1537418130407
}
]
```



**Description** : Fetch all the flows whose destination is the specified IP Address.

**Request URI** : `https://<Ordr_SCE>/Rest/Flows?dstIp=10.200.204.1`

**Output** : [Sample output](#)

**Description** : Fetch all the flows whose source is the specified device.

**Request URI** : `https://192.168.104.182/Rest/Flows?srcMac=00:50:56:07:DB:E4`

**Output** : [Sample output](#)

**Description** : Fetch all the flows whose destination is the specified device.

**Request URI** : `https://<Ordr_SCE>/Rest/Flows?dstMac=52:54:00:01:79:B6`

**Output** : [Sample output](#)

**Description** : Retrieve the current version of the Applications API supported by Ordr SCE.

**Request URI** : `https://<Ordr_SCE>/Rest/Applications/version`

**Output** : Current supported version of /Rest/Applications API is 1.0

**Description** : Fetch all the applications used by the device for the given device MAC.

**Request URI** : `https://<Ordr_SCE>/Rest/Applications?mac=28:63:36:A6:F9:01`

**Output** :

```
{
  "MetaData": {
    "Count": 4
  },
  "Applications": [
    {
      "protocol": "tcp",
      "appName": "ssh",
      "peers": [
        "00:0C:29:3F:A3:20"
      ]
    },
    {
      "protocol": "udp",
      "appName": "mdns",
      "peers": [
        "00:0C:29:02:28:1D",
        "00:0C:29:1C:7B:68",
        "00:0C:29:D7:6D:AB"
      ]
    }
  ]
}
```

```

    ]
  },
  {
    "protocol": "udp",
    "appName": "SNMP",
    "peers": [
      "00:0C:29:02:28:1D",
      "00:0C:29:1C:7B:68",
      "00:0C:29:D7:6D:AB"
    ]
  },
  {
    "protocol": "udp",
    "appName": "ssdp",
    "peers": [
      "9C:93:4E:3C:D7:75",
      "AC:CC:8E:2B:A5:E4"
    ]
  }
]
}

```

**Description** : Fetch all the applications used by the device for the given device IP.

**Request URI** : [https://<Ordr\\_SCE>/Rest/Applications?ip=10.200.205.16](https://<Ordr_SCE>/Rest/Applications?ip=10.200.205.16)

**Output** : [Sample output](#)

**Description** : Fetch all the security incidents detected by the system.

**Request URI** : [https://<Ordr\\_SCE>/Rest/SecurityAlarms](https://<Ordr_SCE>/Rest/SecurityAlarms)

**Output** :

```

{
  "MetaData": {
    "Count": 100,
    "next": "/Rest/SecurityAlarms?clientMacToken=-
2454454372417895144&alarmHashToken=f743d6c2626385ef950754df93824e7df5e10d3e"
  },
  "SecurityAlarms": [
    {
      "alarmHash": "1128e80da32f658cb889938e3e2b5618dd85e568",
      "category": "KNOWN_VULN",
      "categoryType": "KNOWN_VULN",
      "severityLevel": "NORMAL",
      "riskScore": 0,
      "metaData": "",
      "peerId": "NA",
      "recentTimestamp": 1579691452742,
      "incidentType": "FDA-156164:1.5T Signa HDx, 3.0T Signa HDx, 1.5T Signa
HDxt, 3.0T Signa HDxt, Sig",

```

```

    "deviceMac": "00:50:56:D7:97:DB",
    "sensorName": "reports-dpvm-ss48",
    "sensorIp": "172.18.10.16"
  },
  {
    "alarmHash": "a5992a08bb38bc91c8df29107449c24143e93fe6",
    "category": "BAD_URL",
    "categoryType": "URL Malware",
    "severityLevel": "MEDIUM",
    "riskScore": 6,
    "metaData": "http://www.disneylanddaze.com/",
    "peerId": "http://www.disneylanddaze.com/",
    "recentTimestamp": 1579960155038,
    "incidentType": "Malware Site Access",
    "deviceMac": "00:0C:29:89:70:7E",
    "rawVectorGuid": "a-209104--132490175",
    "sensorName": "reports-dpvm-ss48",
    "sensorIp": "172.18.10.16",
    "locationInfo": {
      "country": "United States",
      "countryCode": "US",
      "city": "Bluffdale"
    }
  }
],
<<<<<SNIP>>>>>
{
  "alarmHash": "f86269de68350b1f74cbc8b5df61ea27101fcf35",
  "category": "BAD_IP",
  "categoryType": "Suspicious Traffic",
  "severityLevel": "MEDIUM",
  "riskScore": 6,
  "metaData": "23.129.64.159",
  "peerId": "23.129.64.159",
  "recentTimestamp": 1583330351185,
  "incidentType": "packets to blacklisted destination",
  "deviceMac": "52:54:00:89:82:C5",
  "sensorName": "reports-dpvm-ss48",
  "sensorIp": "172.18.10.16"
}
]
}

```

**Description** : Fetch all the security incidents for the device for the given device MAC.

**Request URI** : <https://192.168.104.182/Rest/SecurityAlarms?mac=<mac>>

**Output** :

```

{
  "MetaData": {
    "Count": 1
  },
  "SecurityAlarms": [
    {

```

```

    "alarmHash": "a5992a08bb38bc91c8df29107449c24143e93fe6",
    "category": "BAD_URL",
    "categoryType": "URL Malware",
    "severityLevel": "MEDIUM",
    "riskScore": 6,
    "metaData": "http://www.disneylanddaze.com/",
    "peerId": "http://www.disneylanddaze.com/",
    "recentTimestamp": 1579960155038,
    "rawVectorGuid": "a-209104--132490175",
    "locationInfo": {
      "country": "United States",
      "countryCode": "US",
      "city": "Bluffdale"
    }
  }
]
}

```

**Description** : Fetch all the security incidents for the device for the given device IP.

**Request URI** : `https://192.168.104.182/Rest/SecurityAlarms?mac=<mac>`

**Output** : [Sample output](#)

**Description** : Fetch all the security incidents for the given alarm category.

**Request URI** : `https://<Ordr_SCE>/Rest/SecurityAlarms?category=BAD_URL`

**Output** :

```

{
  "MetaData": {
    "Count": 1
  },
  "SecurityAlarms": [
    {
      "alarmHash": "a5992a08bb38bc91c8df29107449c24143e93fe6",
      "category": "BAD_URL",
      "categoryType": "URL Malware",
      "severityLevel": "MEDIUM",
      "riskScore": 6,
      "metaData": "http://www.disneylanddaze.com/",
      "peerId": "http://www.disneylanddaze.com/",
      "recentTimestamp": 1579960155038,
      "incidentType": "Malware Site Access",
      "deviceMac": "00:0C:29:89:70:7E",
      "rawVectorGuid": "a-209104--132490175",
      "sensorName": "reports-dpvm-ss48",
      "sensorIp": "172.18.10.16",
      "locationInfo": {
        "country": "United States",
        "countryCode": "US",

```

```

        "city": "Bluffdale"
    }
}
]
}

```

**Description** : Provide a summary of the security alarms by category.

**Request URI** : [https://<Ordr\\_SCE>/Rest/SecurityAlarms/Summary](https://<Ordr_SCE>/Rest/SecurityAlarms/Summary)

**Output** :

```

{
  "summary": [
    {
      "category": "BAD_URL",
      "summary": [
        {
          "categoryType": "URL Malware",
          "severityLevel": "MEDIUM",
          "riskScore": 6,
          "deviceCount": 2
        }
      ]
    },
    {
      "category": "KNOWN_VULN",
      "summary": [
        {
          "categoryType": "KNOWN_VULN",
          "severityLevel": "NORMAL",
          "riskScore": 0,
          "deviceCount": 11
        }
      ]
    },
    {
      "category": "URL_GENERIC",
      "summary": [
        {
          "categoryType": "URL Generic alarm",
          "severityLevel": "MEDIUM",
          "riskScore": 6,
          "deviceCount": 1
        }
      ]
    },
    {
      "category": "DEVICE_SIGNATURE_VIOLATION",
      "summary": [
        {
          "categoryType": "Baseline Flow Violation",
          "severityLevel": "MEDIUM",

```

```
"riskScore": 6,
"deviceCount": 7
}
],
},
{
"category": "PHISHING",
"summary": [
{
"categoryType": "URL Phishing",
"severityLevel": "MEDIUM",
"riskScore": 6,
"deviceCount": 3
}
]
}
]
```

**Description** : Fetch all network equipments.

**Request URI** : [https://<Ordr\\_SCE>/Rest/NetworkDevices](https://<Ordr_SCE>/Rest/NetworkDevices)

**Output** :

```
{
  "networkDeviceInfo": {
    "totalAccessSwitches": 22,
    "totalAccessPorts": 564,
    "totalDot1xDisabledPorts": 564,
    "totalMabDisabledPorts": 564,
    "totalDot1xAndMabDisabledPorts": 564,
    "networkDevices": [
      {
        "name": "HP-3800-48G-PoEP-2SFPP",
        "scrapeIp": "10.100.16.10",
        "accessPorts": 23,
        "dot1xDisabledPorts": 23,
        "mabDisabledPorts": 23,
        "dot1xAndMabDisabledPorts": 23
      },
      {
        "name": "uplink_to_controller",
        "scrapeIp": "10.100.13.3",
        "accessPorts": 23,
        "dot1xDisabledPorts": 23,
        "mabDisabledPorts": 23,
        "dot1xAndMabDisabledPorts": 23
      }
    ]
  }
}
```



# ördr

take control.

[info@ordr.net](mailto:info@ordr.net)  
[www.ordr.net](http://www.ordr.net)



2445 Augustine Drive Suite 601  
Santa Clara, CA 95054

